



Benkenstrasse 254 – CH-4108 Witterswil – www.biolytix.ch

Detection of different *Apis* species with real-time PCR

SOP

Author: M. Schöne-Michling

November 2006

1. ABSTRACT

1.1. PURPOSE

This method describes a practice for detection of desoxynucleinacids (DNA) from different apis species with real-time PCR.

1.2. AREA OF APPLICATION

With the method describe it is possible to detect DNA from *Apis cerana*, *Apis mellifera*, *Apis florea* and *Apis dorsata* with real-time PCR .

1.3. Origin

This method was designed by Mandy Schöne-Michling.

2 PRINCIPLE

The real-time PCR systems are based on detection of mitochondrial rRNA.

3 TEST EQUIPMENT

3.1. Equipment and Material

Important!

Only use analyze clean, for molecular biology grade reagents and water. The handling must be done under sterile conditions.

Powder free gloves have to be used.

- ABI Prism® Sequence Detection System 7700, Applied Biosystems
- Rotor-Gene 3000, Corbett Research
- Mikrocentrifuge
- Pipettes
- Gloves, powder free
- PCR using material

3.2. Reagents

(Use only DNase and Rnase free reagents.)

- Sterile H₂O bidest
- 2x RT PCR Master Mix

| Name of primer and probe | Sequenz 5' – 3' | Target organism |
|--------------------------|-------------------------------------|-----------------------|
| A.cera F | CCTTAGGATAACAGCGTAATAT | <i>Apis cerana</i> |
| A.cera R | CAGACTTAAAATTTTAAACTCCT | |
| A.cera probe | TTGATAGACCACATTGATAAAGATGT (VIC) | |
| A.dors F | CTTAGGATAACAGCGTAATATC | <i>Apis dorsata</i> |
| A.dors R | ATTTAAACTTCTGCATTTAACTTTA | |
| A.dors probe | TTGATAGACCTTATAGATAAAGATGAT (VIC) | |
| A.florea F | GATTAGAAGATATAGAAATGATTT | <i>Apis florea</i> |
| A.florea R | CAATCATCTTTATCAATATGACC | |
| A.florea probe | TAATTTATAAATAAGAATAAATTACCTTA (VIC) | |
| A.meli F | CTTAGGATAACAGCGTAATATC | <i>Apis mellifera</i> |
| A.meli R | AAACTTAAACTACTGCGCCTA | |
| A.meli probe | TAGACCATATAGATAAAGATGTTTGC (FAM) | |

3.3. Stock solutions

all primers: 8 µM primer solution
all probes: 8 µM probe solution

3.4. Reference material

- Reference DNA from *Apis cerana*
- Reference DNA from *Apis dorsata*
- Reference DNA from *Apis florea*
- Reference DNA from *Apis mellifera*

4 EXECUTION

4.1. Preparation of Mastermix solution

| Components | 1 | 10 | Endconcentration |
|----------------------|------|-----|------------------|
| 2x RT PCR Master Mix | 12.5 | 125 | 1 x |
| Forward Primer | 1 | 10 | 320 nM |
| Reverse Primer | 1 | 10 | 320 nM |
| Probe | 0.5 | 5 | 160 nM |
| H ₂ O | 5 | 50 | |
| Summe | 20 | 200 | |

4.2. Adding of DNA extract to the Mastermix solution

- Vortex the mastermix solution after preparation.
- Do 20µl of mastermix solution in a sterile tube for ABI Prism 7700 or Rotor-Gene 3000.
- Add 5µl of template DNA (10 fold solution) to the tubes.
- For a negative control add 5µl of sterile H₂O bidest.
- Close the tubes and spin short.

4.3. Real-time-PCR conditions

| Step | |
|---------------------------|------------------------------------|
| UNG activity | 2 min. / 50 °C |
| Activation of HotStarTaq | 10 min. / 95 °C |
| Amplification (50 cycles) | 15 sec. / 95 °C 60 sec. / 60 °C |

5 ANALYSIS**6 ADVICE**